BIOPATCH™

Protective Disk with CHG
Предпазен диск с CHG
Zaštitni disk sa CHG-om
Ochranný disk s CHG
Beskyttelseskive med CHG
Beschermede schijf met CHG
CHG-ga kaitseketas
Suojalevy ja CHG
Disque de protection avec GCH
Schutzscheibe mit CHG
Προστατευτικός δίσκος με CHG
Védőlemez CHG-vel
Diosca Cosanta le GCH
Disco protettivo con CHG
Aizsargdisks ar CHG
Apsauginis diskas su CHG
Diska Protettiva b'CHG
Dysk ochronny z CHG
Disco protetor com GCH
Disc de protecție cu CHG
Ochranný kotúč s CHG
Zaščitni disk s CHG
Disco protector con CHG
Skyddsskiva med CHG

STERILE EO

ETHICON™
Johnson & Johnson SURGICAL TECHNOLOGIES

CE
2797

UK
CA
0086

Instructions for Use

(Please Read Carefully Before Using)

Device Description

BIOPATCH™ Protective Disk with CHG is a hydrophilic polyurethane absorptive foam with Chlorhexidine Gluconate (CHG). The foam material absorbs up to eight times its own weight in fluid, while the CHG incorporated into the dressing inhibits bacterial growth under the dressing. Chlorhexidine Gluconate is a well-known antiseptic agent with broad-spectrum antimicrobial and antifungal activity.

Performance Characteristics

BIOPATCH™ Protective Disk with CHG is comprised of a white to off-white hydrophilic polyurethane foam impregnated with 200-300 micrograms of chlorhexidine gluconate (CHG) per milligram of foam. The foam is reinforced with a urethane film which is colored blue and printed with text and figures (BIOPATCH™, UP and double arrows). The dressing contains a central aperture with a radial slit to enable placement around a percutaneous device and ensure that the dressing can make 360 degree contact with the skin around the device. Release of CHG from the foam inhibits bacterial growth beneath the dressing and on the skin around the device insertion site. The skin is known to be a major source for microorganisms associated with infectious complications of percutaneous devices. The printed blue film ensures proper placement of the dressing (i.e. blue side up) as the CHG can only release from the side of the dressing without the urethane film. Once placed around the percutaneous device, the dressing is secured to the patient and may remain in place for up to seven days.

Indications for Use/Intended Use

BIOPATCH™ containing Chlorhexidine Gluconate is intended for use as a hydrophilic wound dressing that is used to absorb exudate and to cover a wound caused by the use of vascular and non-vascular percutaneous medical devices such as: IV catheters, central venous lines, arterial catheters, dialysis catheters, peripherally inserted coronary catheters, mid-line catheters, drains, chest tubes, externally placed orthopedic pins, and epidural catheters. It is also intended to reduce local infections, catheter-related blood stream infections (CRBSI), and skin colonization of microorganisms commonly related to CRBSI, in patients with central venous or arterial catheters.

Intended Users

This device is intended to be used by trained health care professionals.

Clinical Benefits

Reduce local infections, catheter-related blood stream infections (CRBSI), and skin colonization of microorganisms commonly related to CRBSI, in patients with central venous or arterial catheters.

Summary of Safety and Clinical Performance (SSCP) is available upon request and on EUDAMED at <https://ec.europa.eu/tools/eudamed>.

Precautions

BIOPATCH™ should not be placed over infected wounds. It is not intended to be used as a treatment of percutaneous device-related infections.

Warnings

WARNING: DO NOT USE BIOPATCH™ ON PREMATURE INFANTS. USE OF THIS PRODUCT ON PREMATURE INFANTS HAS RESULTED IN HYPERSENSITIVITY REACTIONS AND NECROSIS OF THE SKIN.

DO NOT RESTERILIZE/REUSE. REUSE OF THIS DEVICE (OR PORTIONS OF THIS DEVICE) MAY CREATE A RISK OF PRODUCT DEGRADATION AND CROSS-CONTAMINATION, WHICH MAY LEAD TO INFECTION OR TRANSMISSION OF BLOOD-BORNE PATHOGENS TO PATIENTS AND USERS.

FOR EXTERNAL USE ONLY. DO NOT ALLOW THIS PRODUCT TO CONTACT THE EYES, EARS, MOUTH, OR MUCOUS MEMBRANES.

THE SAFETY AND EFFECTIVENESS OF BIOPATCH™ HAS NOT BEEN ESTABLISHED IN CHILDREN UNDER 16 YEARS OF AGE.

DO NOT USE BIOPATCH™ DIRECTLY OVER BURN INJURY OR ON PATIENTS WITH A KNOWN SENSITIVITY TO CHLORHEXIDINE GLUCONATE. ADVERSE REACTIONS TO CHLORHEXIDINE GLUCONATE SUCH AS DERMATITIS, HYPERSENSITIVITY, AND GENERALIZED ALLERGIC REACTIONS ARE VERY RARE, BUT IF ANY SUCH REACTIONS OCCUR, DISCONTINUE USE OF THE DRESSING IMMEDIATELY.

HYPERSENSITIVITY REACTIONS ASSOCIATED WITH THE TOPICAL USE OF CHLORHEXIDINE GLUCONATE HAVE BEEN REPORTED IN SEVERAL COUNTRIES. THE MOST SERIOUS REACTIONS [INCLUDING ANAPHYLAXIS] HAVE OCCURRED IN PATIENTS TREATED WITH LUBRICANTS CONTAINING CHLORHEXIDINE GLUCONATE, WHICH WERE USED DURING URINARY TRACT PROCEDURES. PREPARATIONS OF THIS TYPE ARE NOT APPROVED FOR SALE IN THE U.S. UNDER ANY CIRCUMSTANCES.

CAUTION SHOULD BE USED WHEN USING CHLORHEXIDINE-CONTAINING PREPARATIONS AND THE PATIENT SHOULD BE OBSERVED FOR THE POSSIBILITY OF HYPERSENSITIVITY REACTIONS. THE GOVERNMENT OF JAPAN HAS REPORTED ANAPHYLACTOID-TYPE ADVERSE EVENTS IN 13 PATIENTS WHILE USING CENTRAL VENOUS CATHETERS IMPREGNATED WITH CHLORHEXIDINE.

ANY SERIOUS INCIDENT THAT HAS OCCURRED IN RELATION TO THE DEVICE FOR THE USER AND/OR THE PATIENT SHOULD BE REPORTED TO THE MANUFACTURER AND THE COMPETENT AUTHORITY OF THE MEMBER STATE IN WHICH THE USER AND/OR PATIENT IS ESTABLISHED.

Directions for Use

1. Prepare the skin surrounding the percutaneous device according to hospital protocol.
2. Remove BIOPATCH™ from the sterile package using aseptic technique.
3. Place BIOPATCH™ around the device, making sure the BLUE PRINTED side is facing upward. The WHITE foam side releases the Chlorhexidine Gluconate (CHG) and should be in contact with the patient's skin.
4. In order to ensure easy removal when used with a film dressing, place BIOPATCH™ around the device site in such a way that the device rests upon the slit portion of the BIOPATCH™. The edges of the radial slit must be pushed together and remain in contact to maximize efficacy.
5. Secure the device and BIOPATCH™ to the skin. Ensure complete contact between the skin and BIOPATCH™.
6. Change the patch as necessary, in accordance with facility protocol; dressing changes should occur at a minimum of every 7 days. Dressing changes will be needed more frequently with highly exuding wounds.

7. To remove the transparent film dressing, pick up the corner of the dressing and stretch the dressing away from the device, holding the device in place. (Dressing will partially lift.) Peel back until resistance is felt. Repeatedly stretch and peel as necessary until the dressing is removed.
8. BIOPATCH™ will remain attached to the transparent film dressing, so removal will be simultaneous.

Storage/Transportation Information

- Store between 15°C and 30°C (59°F and 86°F).
- It is to be stored in its original packaging.
- Expiration date of the product is indicated as year (4 digits), month (2 digits) and day (2 digits).
- Do not resterilize. Discard all open and unused portions of the device.
- Do not use if the package is opened or damaged. Do not use if seal is broken or compromised.
- After use, handle and dispose of all unused product and packaging in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

NOTE: Over time, the BIOPATCH™ may turn yellow in color. This coloration does not reduce the antimicrobial efficacy of the dressing.

How Supplied

BIOPATCH™ is supplied sterile. Each package contains a single disk. BIOPATCH™ is intended for single use only.

Do not resterilize.

| Product Code Ending | BIOPATCH™ Size | Maximum Amount of CHG per Dressing |
|---------------------|--|------------------------------------|
| -150 | 1" DISK (2.5 cm) 4.0 mm center hole | 92 mg |
| -151 | 3/4" DISK (1.9 cm) 1.5 mm center hole | 53 mg |
| -152 | 1" DISK (2.5 cm) 7.0 mm center hole | 86.8 mg |

Product Disposal

Follow approved medical waste disposal methods as specified by your patient care facility or local regulations.

Incident Reporting

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the local competent authority of the country in which the user and/or patient is established.

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Integra LifeSciences Corporation
Ethicon Inc.
ILS Services Switzerland Ltd
Integra NeuroSciences Ltd