



Health Care
Doctor
Hospital
Pharmacist
Nurse
Dentist
First Aid
Surgeon
Emergency

HEALTHCARE, INC.

MANUFACTURER OF QUALITY
MEDICAL DEVICES AND
DISPOSABLE KITS

DIAGNOSTIC PROCEDURE TRAYS

***Where Quality Exceeds
Expectations!***

8(a) Certified Small Minority Business Manufacturer, by U.S.
Small Business Administration

**Niles, Illinois 60714
Tel: (888) CALL-LSL**

**Visit Our Website @
www.lslhealthcare.com**



TABLE OF CONTENTS

Lumbar Puncture Trays	3-4
Myelogram	5
Arthrogram	6
Amniocentesis Trays	7
Paracentesis / Thoracentesis Trays	8-9
Bone Marrow aspiration & Biopsy Trays	10-11
Soft Tissue Biopsy Trays	12
Basic Diagnostic Trays	13





LUMBAR PUNCTURE TRAYS

LSL's lumbar puncture trays contain the core components needed for any lumbar puncture procedure. Our lumbar puncture tray offerings include a traditional and safety needle version, with or without a manometer as well as an option for the adult or pediatric patient.

LSL always goes a step further by allowing you to customize a kit specific to your needs. Add in extra gauze, switch out a smaller bandage for a larger bandage; whatever your requirements, you know you can trust our kits to contain the highest quality products from the industry's leading manufacturers.

CORE COMPONENTS

COMPONENTS

QTY

Adhesive bandage 1 x 1 1/2 in.	1
Drape, fenestrated with tape	1
Extension tubing, 5 in	1
Foam needle block	1
Gauze, 8-ply, 4 x 4 in	10
Lidocaine 1% ampoule, 5 mL	1
Medication label sheet.	1
PVP-I swabstick, 3/pk	1
Stopcock, 3-way	1
Syringe, luer-lock, 5 mL	1
Towel, 17 x 19 in.	1
Utility marker	1
Vial with cap, clear, #1, 10 mL	1
Vial with cap, clear, #2, 10 mL	1
Vial with cap, clear, #3, 10 mL	1
Vial with cap, clear, #4, 10 mL	1
Wrap, 30 x 30 in.	1

ACCESSORY SET

VIAL SET, 10 mL COMPONENTS

QTY

Vial with cap, clear, #1, 10	1
Vial with cap, clear, #2, 10	1
Vial with cap, clear, #3, 10	1
Vial with cap, clear, #4, 10	1



CONVENTIONAL KITS

**Lumbar puncture w/ spinal needle,
20g x 3 ½ in.**

Reorder #: 6063

QTY

Filter Needle, 19g x 1 ½ in.
Manometer, two-piece
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Spinal Needle, 20g x 3 ½ in.
Lumbar puncture w/ spinal needle,
pediatric, 22g x 1 ½ in.

1
1
1
1
1

Reorder #: 6083

QTY

Filter Needle, 19g x 1 ½ in.
Manometer, one-piece
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Spinal Needle, 22g x 1 ½ in.
Lumbar puncture w/ spinal needle,
22g x 3 ½ in.

1
1
1
1
1

Reorder #: 6088

QTY

Filter Needle, 19g x 1 ½ in.
Manometer, two-piece
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Spinal Needle, 22g x 3 ½ in.
Lumbar puncture atraumatic w/
spinal needle, 22g x 3 ½ in.

1
1
1
1
1

Reorder #: 6093

QTY

Filter Needle, 19g x 1 ½ in.
Manometer, two-piece
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Spinal Needle, 22g x 3 ½ in.
Lumbar puncture, no manometer

1
1
1
1
1

Reorder #: 6096

QTY

Filter Needle, 19g x 1 ½ in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Spinal Needle, 20g x 3 ½ in.
Lumbar puncture w/ spinal needle,
18g x 3 ½ in.

1
1
1
1
1

Reorder #: 6098

QTY

Filter Needle, 19g x 1 ½ in.
Manometer, two-piece
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Spinal Needle, 18g x 3 ½ in.

1
1
1
1
1

SAFETY KITS

**Safety Lumbar puncture w/ spinal needle,
20g x 3 ½ in.**

Reorder #: 6063SF

QTY

Filter Straw, 5 Micron, 1 ¾ in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Manometer, two-piece
Spinal Needle, 20g x 3 ½ in.
Safety Lumbar puncture w/ spinal
needle, pediatric, 22g x 1 ½ in.

1
1
1
1
1

Reorder #: 6083SF

QTY

Filter Straw, 5 Micron, 1 ¾ in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Manometer, one-piece
Spinal Needle, 22g x 1 ½ in.
Safety Lumbar puncture w/ spinal
needle, 22g x 3 ½ in.

1
1
1
1
1

Reorder #: 6088SF

QTY

Filter Straw, 5 Micron, 1 ¾ in.
Manometer, two-piece
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Spinal Needle, 22g x 3 ½ in.
Safety Lumbar puncture atraumatic
w/spinal needle, 22g x 3 ½ in.

1
1
1
1
1

Reorder #: 6093SF

QTY

Filter Straw, 5 Micron, 1 ¾ in.
Manometer, two-piece
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Spinal Needle, 22g x 3 ½ in.
Safety Lumbar puncture w/ spinal
needle, 18g x 3 ½ in.

1
1
1
1
1

Reorder #: 6098SF

QTY

Filter Needle, 19g x 1 ½ in.
Manometer, two-piece
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Spinal Needle, 18g x 3 ½ in.

1
1
1
1
1

All products listed on this page are
not made with natural rubber latex.
Needles are manufactured by either
BD or Covidien





MYELOGRAM TRAYS

LSL's myelogram trays contain the fundamental components needed to perform a myelogram procedure. We offer a conventional or safety version and customization is always an option with LSL trays.

CORE COMPONENTS

COMPONENTS

QTY

Adhesive bandage 1 x 1 ½ in.	1
Drape, fenestrated with tape	1
Extension tubing, 5 in	1
Foam needle block	1
Gauze, 8-ply, 4 x 4 in	10
Lidocaine 1% ampoule, 5 mL	1
Medication label sheet	1
PVP-I swabstick, 3/pk	1
Stopcock, 3-way	1
Syringe, luer-lock, 5 mL	1
Towel, 17 x 19 in.	1
Utility marker	1
Vial with cap, clear, #1, 10 mL	1
Vial with cap, clear, #2, 10 mL	1
Wrap, 30 x 30 in.	1

CONVENTIONAL KITS

Myelogram with spinal needle, 22g x 3 1/2 in.

Reorder #: 6026

QTY

Filter Needle, 19g x 1 ½ in.	1
Needle, 22g x 1 ½ in.	1
Needle, 25g x 5/8 in.	1
Spinal Needle, 22g x 3 ½ in.	1
Myelogram, atraumatic with spinal needle, 22g x 3 1/2 in.	

Reorder #: 6028

QTY

Filter Needle, 19g x 1 ½ in.	1
Needle, 22g x 1 ½ in.	1
Needle, 25g x 5/8 in.	1
Spinal Needle, 22g x 3 ½ in.	1

SAFETY KITS

Safety myelogram with spinal needle, 22g x 3 1/2 in.

Reorder #: 6026SF

QTY

Filter Straw, 5 Micron, 1 ¾ in.	1
Needle, 22g x 1 ½ in.	1
Needle, 25g x 5/8 in.	1
Spinal Needle, 22g x 3 ½ in.	1

All products listed on this page are not made with natural rubber latex.
Needles are manufactured by either BD or Covidien





ARTHROGRAM TRAYS

LSL's arthrogram trays include all of the core essentials necessary for arthrography procedures. LSL partners with only top manufacturers to ensure our kits contain only the highest quality products. Conventional and safety versions are available.

CORE COMPONENTS

COMPONENTS

QTY

Adhesive bandage 1 x 1 ½ in.	1
ChloraPrep® applicator, 3.0 mL	1
Drape, fenestrated with tape	1
Extension tubing, 20 in	1
Foam needle block	1
Gauze, 8-ply, 4 x 4 in	10
Lidocaine 1% ampoule, 5 mL	1
Medication label sheet	1
Fill Needle, 18g x 1 1/2 in.	1
Foam needle block	1
Syringe, luer-lock, 5 mL	1
Syringe, luer-lock, 20 mL	1
Towel, 17 x 19 in.	1
Utility marker	1
Wrap, 30 x 30 in.	1

CONVENTIONAL KITS

Arthrogram with spinal needle,

Reorder #: 6038

QTY

Filter Needle, 19g x 1 ½ in.	1
Needle, 22g x 1 ½ in.	1
Needle, 25g x 5/8 in	1

SAFETY KITS

Safety Arthrogram with spinal needle,

Reorder #: 6038SF

QTY

FILTER STRAW®, 5	1
Micron, 1 ¾ in.	
Needle, 22g x 1 ½ in.	1
Needle, 25g x 5/8 in.	1

All products listed on this page are not made with natural rubber latex.
Needles are manufactured by either BD or Covidien





AMNIOCENTESIS TRAYS

LSL offers amniocentesis trays that feature advanced, innovative components. Safety is paramount for Mom and Baby, that is why our kits include the most innovative products including the echogenic needle, allowing the clinician to view the needle under ultrasound. Tailor your amniocentesis tray to meet your needs. Conventional and Safety versions available.

CORE COMPONENTS

COMPONENTS

QTY

Adhesive bandage 1 x 1 ½ in.	1
Cap, male/female, luer-lock adapter	1
ChloraPrep® applicator, 3.0 mL	1
Drape, fenestrated with tape	1
Echogenic procedure needle, 20 g x 3 ½ in.	1
Extension tubing, 5 in	1
Gauze, 8-ply, 4 x 4 in	10
Lidocaine 1% ampoule, 5 mL	1
Medication label sheet	1
Foam needle block	1
Ruler 6 in.	1
Syringe, luer-lock, 5 mL	1
Syringe, luer-lock, 20 mL	1
Towel, 17 x 19 in.	1
Utility marker	1
Vial with cap, amber, #1, 10 mL	1
Vial with cap, amber, #2, 10 mL	1
Vial with cap, amber, #3, 10 mL	1
Wrap, 30 x 30 in.	1

CONVENTIONAL KITS

**Amniocentesis with echogenic needle,
20g x 3 ½ in.**

Reorder #: 6057

QTY

Filter Needle, 19g x 1 ½ in.	1
Needle, 22g x 1 ½ in.	1
Needle, 25g x 5/8 in.	1

SAFETY KITS

**Safety amniocentesis with echogenic
needle, 20g x 3 1/2 in.**

Reorder #: 6057SF

QTY

FILTER STRAW®, 5 Micron, 1 ¾ in.	1
Needle, 22g x 1 ½ in.	1
Needle, 25g x 5/8 in.	1

All products listed on this page are not made with natural rubber latex.
Needles are manufactured by either BD or Covidien





PARACENTESIS / THORACENTESIS TRAYS

LSL's paracentesis/thoracentesis trays contain cutting-edge technology allowing for the rapid removal of fluid buildup from the pleural or abdominal cavity. LSL kits are customizable and offer conventional or safety option.

CORE COMPONENTS

COMPONENTS

QTY

#11—Safety Scalpel	1
Adhesive bandage 1 x 1 ½ in.	1
Cap, male/female, luer-lock adapter	1
ChloraPrep® applicator, 3.0 mL	1
Drape, fenestrated with tape	1
Extension tubing with roller clamp and 16 g	1
Gauze, 8-ply, 4 x 4 in	10
Lidocaine 1% ampoule, 5 mL	1
Medication label sheet	1
Foam needle block	1
Syringe, luer-lock, 5 mL	1
Syringe, luer-lock, 30 mL	1
Towel, 17 x 19 in.	1
Transparent dressing, 4 x 4 ¾ in.	1
Utility marker	1
Wrap, 30 x 30 in.	1

ACCESSORY SET

Centesis drainage bag set

Reorder #: 6017

QTY

Drainage Bag	1
Drainage Bag, Tubing Set	1

All products listed on this page are not made with natural rubber latex.
Needles are manufactured by either BD or Covidien





PARACENTESIS / THORACENTESIS TRAYS

CONVENTIONAL KITS

Centesis, 8 french x 5 in.

Reorder #: 6019

QTY

Filter Needle, 19g x 1 ½ in.
Fluid drainage system 8 french
x 5 in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.

Centesis with aspirating needle,
17g x 6 in.

1
1
1
1
1

Reorder #: 6022

QTY

Extension tubing with stopcock, 3 in.
Filter Needle, 19g x 1 ½ in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Thoracentesis needle, 14g Catheter
over 17 g x 6 in. needle

Centesis, no catheter / needle

1
1
1
1
1

Reorder #: 6023

QTY

Adapter
Drainage bag
Drainage bag Tubing Set
Filter Needle, 19g x 1 ½ in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Spinal Needle, 22g x 3 ½ in.

1
1
1
1
1
1
1

SAFETY KITS

Safety Centesis, 8 french x 5 in.

Reorder #: 6019SF

QTY

FILTER STRAW®, 5 Micron,
1 ¾ in.
Safety Fluid drainage system
8 french x 5 in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Centesis with aspirating needle,
17g x 6 in.

1
1
1
1

Reorder #: 6022SF

QTY

FILTER STRAW®, 5 Micron,
1 ¾ in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Thoracentesis needle, 14g Catheter
over 17 g x 6 in. needle
Safety Centesis, no catheter /
needle

1
1
1
1

Reorder #: 6023SF

QTY

Adapter
Drainage bag
Drainage bag Tubing Set
FILTER STRAW®, 5 Micron,
1 ¾ in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.

1
1
1
1
1
1

All products listed on this page are not made with natural rubber latex.
Needles are manufactured by either BD or Covidien





BONE MARROW ASPIRATION & BIOPSY TRAYS

LSL's bone marrow aspiration and biopsy trays provide the options needed to collect samples from various sample sites. You can customize our tray to include a J-type or I-type (or both) bone marrow aspiration needle in either conventional or safety versions.

CORE COMPONENTS

COMPONENTS

QTY

#11—Safety Scalpel	1
Adhesive bandage 1 x 1 ½ in.	1
Chloraprep® applicator, 3.0 mL	1
Drape, fenestrated with tape	1
Foam needle block	1
Gauze, 8-ply, 4 x 4 in	10
Lidocaine 1% ampoule, 5 mL	1
Medication label sheet	1
Syringe, luer-lock, 10 mL	1
Syringe, luer-lock, 20 mL	1
Towel, 17 x 19 in.	1
Utility marker	1
Wrap, 30 x 30 in.	1

SINGLE NEEDLE OPTIONS

Needle type

Reorder #:

Type

QTY

7465	J-type	8g x 4 in.
7467	J-type	11g x 4 in.
7458	J-type	11g x 6 in.
7469	J-type	13g x 2 ½ in.
7473	I-type	15g

All products listed on this page are not made with natural rubber latex.
Needles are manufactured by either BD or Covidien



CONVENTIONAL KITS

Bone marrow biopsy with J-type needle, 11g x 4 in.

Reorder #: 7483

QTY

Bone marrow biopsy extraction cannula
Bone marrow biopsy probe
Bone marrow biopsy probe guide
Bone marrow biopsy, J-Type, 11g x 4 in.
Filter Needle, 19g x 1 ½ in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Bone marrow biopsy with J-type and I-type needles

1
1
1
1
1
1
1

Reorder #: 7484

QTY

Bone marrow biopsy extraction cannula
Bone marrow biopsy probe
Bone marrow biopsy probe guide
Bone marrow biopsy, I-type, with adjustable guard, 15g
Bone marrow biopsy, J-Type, 11g x 4 in.
Filter Needle, 19g x 1 ½ in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Bone marrow biopsy with I-type needle, 15g

1
1
1
1
1
1
1
1

Reorder #: 7485

QTY

Bone marrow biopsy, I-type, with adjustable guard, 15g
Filter Needle, 19g x 1 ½ in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Bone marrow biopsy, no needle

1
1
1
1

Reorder #: 7486

QTY

Filter Needle, 19g x 1 ½ in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Bone marrow biopsy with J-type and I-type needles

1
1
1

Reorder #: 7487

QTY

Bone marrow biopsy extraction cannula
Bone marrow biopsy probe
Bone marrow biopsy probe guide
Bone marrow biopsy, I-type, with adjustable guard, 15g
Bone marrow biopsy, J-Type, 8g x 4 in.
Filter Needle, 19g x 1 ½ in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.

1
1
1
1
1
1
1
1

SAFETY KITS

Safety bone marrow biopsy with J-type needle, 11g x 4 in.

Reorder #: 6063SF

QTY

Bone marrow biopsy extraction cannula
Bone marrow biopsy probe
Bone marrow biopsy probe guide
Bone marrow biopsy, J-Type, 11g x 4 in.
FILTER STRAW®, 5 Micron, 1 3/4 in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Safety bone marrow biopsy, J-type and I-type needles

1
1
1
1
1
1
1

Reorder #: 7484SF

QTY

Bone marrow biopsy extraction cannula
Bone marrow biopsy probe
Bone marrow biopsy probe guide
Bone marrow biopsy, I-type, with adjustable guard, 15g
Bone marrow biopsy, J-Type, 11g x 4 in.
FILTER STRAW®, 5 Micron, 1 ¾ in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Safety bone marrow biopsy, no needle

1
1
1
1
1
1
1
1

Reorder #: 7486SF

QTY

FILTER STRAW®, 5 Micron, 1 ¾ in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.
Safety bone marrow biopsy with J-type and I-type needles

1
1
1

Reorder #: 7487SF

QTY

Bone marrow biopsy extraction cannula
Bone marrow biopsy probe
Bone marrow biopsy probe guide
Bone marrow biopsy, I-type, with adjustable guard, 15g
Bone marrow biopsy, J-Type, 8g x 4 in.
FILTER STRAW®, 5 Micron, 1 3/4 in.
Needle, 22g x 1 ½ in.
Needle, 25g x 5/8 in.

1
1
1
1
1
1
1
1

All products listed on this page are not made with natural rubber latex. Needles are manufactured by either BD or Covidien





SOFT TISSUE BIOPSY TRAYS

LSL's soft tissue biopsy trays contain the key essentials needed to handle your soft tissue biopsy procedures. Available in conventional or safety versions.

CORE COMPONENTS

COMPONENTS

QTY

#11—Safety Scalpel	1
Adhesive bandage 1 x 1 ½ in.	1
Blunt Fill Needle, 18 g x 1 ½ in.	1
Cap, male/female, luer-lock adapter	1
ChloraPrep® applicator, 3.0 mL	1
Drape, fenestrated with tape	1
Foam needle block	1
Gauze, 8-ply, 4 x 4 in	10
Lidocaine 1% ampoule, 5 mL	1
Medication label sheet	1
Ruler 6 in.	1
Syringe, luer-lock, 3 mL	1
Syringe, luer-lock, 5 mL	1
Syringe, luer-lock, 20 mL	1
Towel, 17 x 19 in.	1
Utility marker	1
Vial with cap, clear, #1, 10 mL	1
Vial with cap, clear, #2, 10 mL	1
Wrap, 30 x 30 in.	1

CONVENTIONAL KITS

Soft tissue biopsy

Reorder #: 7488

QTY

Filter Needle, 19g x 1 ½ in.	1
Needle, 22g x 1 ½ in.	1
Needle, 25g x 5/8 in.	1

SAFETY KITS

Safety soft tissue biopsy

Reorder #: 7488SF

QTY

FILTER STRAW®, 5	1
Micron, 1 ¾ in.	
Needle, 22g x 1 ½ in.	1
Needle, 25g x 5/8 in.	1

All products listed on this page are not made with natural rubber latex.
Needles are manufactured by either BD or Covidien





BASIC DIAGNOSTIC TRAYS

Let's start with the basics. LSL's basic diagnostic trays include all the necessary components to handle a wide range of diagnostic procedures. Customize your tray with best-in-class products that meet today's safety requirements.

CORE COMPONENTS

COMPONENTS

QTY

#11—Safety Scalpel	1
Adhesive bandage 1 x 1 ½ in.	1
Chloraprep® applicator, 3.0 mL	1
Drape, fenestrated with tape	1
Foam needle block	1
Gauze, 8-ply, 4 x 4 in	10
Lidocaine 1% ampoule, 5 mL	1
Medication label sheet	1
Syringe, luer-lock, 5 mL	1
Towel, 17 x 19 in.	1
Utility marker	1
Wrap, 30 x 30 in.	1

CONVENTIONAL KITS

Basic diagnostic tray

Reorder #: 6006

QTY

Filter Needle, 19g x 1 ½ in.	1
Needle, 22g x 1 ½ in.	1
Needle, 25g x 5/8 in.	1

SAFETY KITS

Safety diagnostic tray

Reorder #: 6006SF

QTY

FILTER STRAW®, 5 Micron, 1 ¾ in.	1
Needle, 22g x 1 ½ in.	1
Needle, 25g x 5/8 in.	1

All products listed on this page are not made with natural rubber latex.
Needles are manufactured by either BD or Covidien

